## Dental Instruments Instructions For Use

### Indications

Intra-Lock® instruments are used for dental implant procedures such as site development, implant placement, and implant restorations within the specific indications of each implant system. The label on each instrument contains important product information including whether the instrument is supplied sterile or non-sterile and if the instrument is for single-use only. Instruments and kits supplied non-sterile must be cleaned and sterilized prior to first use and each use thereafter.

Intra-Lock® instruments and surgical kits are indicated for use in site development, placement, and restoration of Intra-Lock® implants and associated components. The Intra-Lock® Single Use Drills are indicated for use with the Intra-Lock® MDL® implants only. No attempt should be made to re-use or re-sterilize the drills. Intra-Lock® assumes no responsibility for attempted re-use or re-sterilization between patients.

### CONTRAINDICATIONS

Intra-Lock® instruments should not be used with patients that have allergies to the specific materials used, including stainless steel and titanium alloy.

### Material

Dental instruments are manufactured from Stainless Steel and Titanium 6AL-4V Alloy.

### How Supplied

#### Sterile Dental Drills

Dental drills that are provided sterile (by gamma radiation) and are intended for single patient use only.

The label on the package provides the lot number, product description, catalog reference number and expiration date. To ensure sterility, dental drills must be used before the end of the expiration date indicated on the outer package label.

Prior to using dental drills, inspect the package and labeling for integrity. If the package has been opened or damaged or the device has been contaminated or exceeded the expiration date before the first use, return dental drill for replacement.

#### Non-Sterile Dental Drills and Instruments

Non-sterile drills MUST be steam sterilized prior to use. See Steam sterilization instructions. The label on the package provides the lot number, product description and catalog reference number.

Prior to using dental drills, inspect the package and labeling for integrity. If the new package has been opened or damaged follow the instructions for steam sterilization.

## Directions for Use

Proper surgical procedures and restorative techniques are the responsibility of the medical professional. Each clinician must evaluate the appropriateness of the procedure used based on personal medical training and experience as applied to the patient case at hand. Intra-Lock® strongly recommends completion of dental implant courses and strict adherence to the instructions pertaining to Intra-Lock® products.

Drills and Taps: Should be replaced when wear is noticed, such as a decrease in cutting efficiency or when signs of discoloration appear. Cutting instruments should be replaced after approximately 12 to 20 osteotomy cycles, depending on bone density. Intra-Lock® recommends use of a drill usage chart to track drill use and to ensure drills are replaced as directed. Single Use Drills: Insert the drill into a handpiece and twist protective tubing off prior to use.

## Handling and Steam Sterilization

Always handle the product with powder-free gloves and avoid contact with hard objects that may damage the surface. If the product is supplied sterile, it should be considered sterile unless the package has been opened or damaged. Using accepted sterile technique, remove product from the package only after the correct size has been determined. If the product is supplied non-sterile or when reprocessing: Remove and discard any shipping material before initial sterilization. Disassemble the surgical tray before each cleaning. Individual instruments must be removed from packaging and cleaned prior to sterilization either within a tray or individually. All instruments must be thoroughly cleaned and sterilized. The following cleaning protocol must be used:

- 1) Prepare a detergent bath in a sterile container using a broad-spectrum cleaning or disinfecting agent such as Hu-Friedy's Enzymax® per the manufacturer's recommendations.
- 2) Brush the product to remove visible debris using a soft bristled brush moistened with the prepared detergent solution.
- 3) Thoroughly rinse product under running utility tap water.
- 4) Place product in the sterile container filled with the prepared detergent solution and sonicate for ten (10) minutes.
- 5) Thoroughly rinse product under running utility tap water.
- 6) Spray or wipe product with 70% IPA.
- 7) Blot product dry with clean lint free cloth.

For surgical kits, place the instruments in the appropriate locations in the surgical tray. For sterilization, place the cleaned kit in an approved sterilization bag or wrap and run through one of the following qualified sterilization cycles:

Sterilization Cycles			
<b>Reference:</b>	ANSI/AAMI TIR12:2010	ANSI/AAMI TIR12:2010	UK HTM 01-01 Part C:2016
<b>Type:</b>	Gravity Steam	Prevacuum Steam	Prevacuum Steam
<b>Exposure Time and Temperature:</b>	15 minutes at 132°C (270°F)	4 minutes at 132°C (270°F)	3 minutes at 134°C (273°F)
<b>Minimum Dry Time:</b>	30 minutes	20 minutes	20 minutes

For individual instruments, place the cleaned instrument in an approved sterilization bag or wrap and run through the following qualified sterilization cycle:

Sterilization Cycles	
<b>Reference:</b>	ANSI/AAMI TIR12:2010
<b>Type:</b>	Gravity Steam
<b>Exposure Time and Temperature:</b>	30 minutes at 121°C (250°F)
<b>Minimum Dry Time:</b>	30 minutes

It is recommended to include a 30-minute cool-down period before removing the product from the sterilization bag or wrap.

 **Warning**

Improper cleaning may lead to inadequate sterilization. Failure to completely dry instruments during autoclaving may leave moisture and cause discoloration and oxidation. The use of hydrogen peroxide or other oxidizing agents will damage the surface of the instruments. Periodic testing, cleaning, and calibration of the autoclave equipment is recommended to ensure the unit remains in proper working order.

## Warnings

- Steam Sterilization instructions have not been validated to reprocess dental drills for use with multiple patients. After use on any patient, the dental drill may include extensive tissue and/or bone debris and therefore require more extensive cleaning and sterilization techniques.
- If dental drills are reused, the user or reprocessing facility must validate the parameters they utilize, which must include an inspection for potential wear and damage incurred to the dental drill from the previous use.

## Drilling Technique

### Warnings

- Osseous drilling can generate detrimental heat that may cause thermal bone necrosis and subsequent failure of the dental implant to osseointegrate. It is critically important that all bone-cutting procedures are CONDUCTED AT SLOW SPEEDS. Profuse, internal, and/or external irrigation is mandatory.
- The slow-speed, highly irrigated drilling procedure is conducted while angling the dental drill such that the direction of the drill bisects the ridge. The dental drill should also be held vertically, avoiding a mesial or distal cant.
- Depth gauge/alignment components can be periodically inserted into the osteotomy site to monitor the angle of penetration. Successively larger drills are used until the desired diameter is achieved.
- Drilling is performed with a precise, up and down pumping action. The drill angle is maintained in order to preserve the concentricity of the hole, while the pumping action allows for incremental depth penetration and periodic cleansing of the flutes.
- It is recommended to drill a Pilot Hole when the bone is of high density using a “Lancer” Pilot drill extending to a depth of 8 mm.

 **Warning**

***Caution: Federal law restricts this device to sale by or on the order of a licensed dentist or physician.***

 **Warning**

***Improper cleaning may lead to inadequate sterilization. Failure to completely dry instruments during autoclaving may leave moisture and cause discoloration and oxidation. The use of hydrogen peroxide or other oxidizing agents will damage the surface of the instruments. Periodic testing, cleaning, and calibration of the autoclave equipment is recommended to ensure the unit remains in proper working order.***

Intra-Lock® International Inc.  
6560 West Rogers Circle, Bldg. 24  
Boca Raton FL 33487 – USA  
[www.intra-lock.com](http://www.intra-lock.com)

Intra-Lock System Europa, S.p.A.  
Via Fabrizio Pinto 16  
Salerno, Salerno, 84124 Italy  
+ 39 089 233 045

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