## Dental Prosthetics Instructions For Use

### Indications

Intra-Lock® abutments are intended for use with Intra-Lock® dental implants to support a prosthetic device in partially or fully edentulous patients. The abutments may be used in single and/or multiple tooth application in the mandible or maxilla.

Intra-Lock® UCLA Gold Abutment Cylinders are used with Intra-Lock® dental implants to create custom prosthetics in partially or fully edentulous patients and can be used in the mandible or maxilla. The Abutments may be used in single and/or multiple tooth application in the mandible or maxilla.

### Contraindications

Patients with alcohol addiction or psychiatric disorders, blood dyscrasias, uncontrolled diabetes, hyperthyroidism, oral infections, malignancies, or patients who have had myocardial infarction within the last 12 months.

Patients with systemic diseases that compromise the immune system, such as AIDS, patients on medications that would compromise healing of an implant site, patients with a history of poor or noncompliance to oral hygiene procedures, or patients who cannot maintain oral hygiene procedures if implants are placed.

Tobacco usage, bruxism, nail biting, and pencil biting increases the occurrence of complications and failures.

Additionally, do not use with patients who have a hypersensitivity to one or more of the metals contained in the alloy.

### Material

Intra-Lock® abutments are manufactured from Titanium 6AL-4V ELI Alloy, ASTM F136 or Commercially Pure Titanium, ASTM F-67. Intra-Lock® manufactures a temporary abutment from PEEK and Titanium 6AL-4V ELI.

Intra-Lock® UCLA Gold Abutment Cylinders are from high gold alloy in the following composition percentages: 58% Au (gold), 22% Pt (platinum), 19% Pd (palladium), 1% Ir (iron) Technical data is available for gold alloy material upon request.

The dental implants and retention screws are manufactured from titanium alloy (Ti 6Al-4V) to ASTM F136.

### How Supplied

Intra-Lock® dental abutments are provided non-sterile and intended for single use only. Abutments should be sterilized prior to use by the instructions provided. Prior to using Intra-

Lock® dental abutments, inspect the package and labeling for integrity. If the device is opened, damaged, or contaminated in any way, it must not be used.

Intra---Lock® UCLA Gold Abutment Cylinders are non-sterile and intended for single use only. Abutments should be sterilized prior to use by the instructions provided. Prior to using Intra---Lock® UCLA Gold Abutment Cylinders, inspect packaging and labeling for integrity. If the device is opened, damaged, or contaminated in any way, DO NOT USE.

### Handling and Sterilization

Always handle the product with powder-free gloves and avoid contact with hard objects that may damage the surface. Remove product from the package only after the correct size has been determined. Non-sterile product must be cleaned and sterilized. The following cleaning protocol must be used:

- 1) Prepare a detergent bath in a sterile container using a broad-spectrum cleaning or disinfecting agent such as Hu-Friedy's Enzymax® per the manufacturer's recommendations.
- 2) Brush the product to remove visible debris using a soft bristled brush moistened with the prepared detergent solution.
- 3) Thoroughly rinse product under running utility tap water.
- 4) Place product in the sterile container filled with the prepared detergent solution and sonicate for two (2) minutes.
- 5) Thoroughly rinse product under running utility tap water.
- 6) Spray or wipe product with 70% IPA.
- 7) Blot product dry with clean lint free cloth.

For sterilization of non-sterile product (not including ceramic prosthetics), place product in an approved sterilization bag or wrap and run through one of the following qualified sterilization cycles:

<b>Sterilization Cycles</b>				
<b>Reference:</b>	ANSI/AAMI TIR12:2010	ANSI/AAMI TIR12:2010	ANSI/AAMI TIR12:2010	UK HTM 01-01 Part C:2016
<b>Type:</b>	Gravity Steam	Gravity Steam	Prevacuum Steam	Prevacuum Steam
<b>Exposure Time and Temperature:</b>	30 minutes at 121°C (250°F)	15 minutes at 132°C (270°F)	4 minutes at 132°C (270°F)	3 minutes at 134°C (273°F)
<b>Minimum Dry Time:</b>	30 minutes	30 minutes	20 minutes	20 minutes

For non-sterile ceramic prosthetics, dry heat sterilize at 170°C (338°F) for sixty (60) minutes.

It is recommended to include a 30-minute cool-down period before removing the product from the sterilization bag or wrap.

### Processing Information for Gold Abutment Cylinder

Laboratory protocol: Fabricate custom abutments by creating a wax-up over the Intra-Lock® UCLA Gold Abutment Cylinders in preparation for a standard investing technique. Maintain the cylinder at the following stabilization temperature ranges.

Conventional technique alloys: between 650 & 700°C. (1200 & 1292°F)

Ceramic alloys: between 800 & 850°C (1472 & 1562°F) for high gold alloys, between 850 & 900°C (1562 & 1652°F) for Palladium based Precious alloys.

Proceed with casting, bringing the alloy to a temperature 150°C over the liquids.

### Warnings

- Dental implant surgery is a complex dental procedure. Appropriate and adequate training in all phases of implant procedures and proper technique is strongly recommended prior to implant use.
- Improper patient selection, diagnosis, treatment planning or technique can result in implant failure and/or loss of supportive bone.
- Care must be taken if performing electrosurgery around a dental implant. Electrosurgery generates heat, which can be conducted through a metallic implant and cause damage to surrounding tissue or bone.
- The external surface of Intra-Lock® Dental Implants should only come in contact with titanium surfaced instruments.


### Warning

**Caution: Federal law restricts this device to sale by or on the order of a licensed dentist or physician.**

Intra-Lock® International Inc.  
6560 West Rogers Circle, Bldg. 24  
Boca Raton FL 33487 – USA  
[www.intra-lock.com](http://www.intra-lock.com)

Intra-Lock System Europa, S.p.A.  
Via Fabrizio Pinto 16  
Salerno, Salerno, 84124 Italy  
+ 39 089 233 045

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