|                  |                                                 |    |
|------------------|-------------------------------------------------|----|
| <b>English</b>   | Instructions for Use: Dental Prosthetics        | 2  |
| <b>Español</b>   | Instrucciones de uso: Prótesis Dentales         | 7  |
| <b>русский</b>   | Инструкция по применению: Зубные протезы        | 12 |
| <b>Deutsch</b>   | Gebrauchsanweisung: Dentale Prothetik           | 18 |
| <b>Français</b>  | Instructions d'utilisation: Prothèses Dentaires | 24 |
| <b>Italiano</b>  | Istruzioni per l'uso: Protesi Dentali           | 29 |
| <b>Português</b> | Instruções de utilização: Próteses Dentárias    | 34 |
| <b>Türkçe</b>    | Kullanım talimatları: Dental Prostatik Ürünler  | 39 |
| <b>简体中文</b>      | 使用说明：牙科义齿                                       | 44 |
| <b>日本語</b>       | 使用説明書：歯科プロテーゼ                                   | 49 |
| <b>한국어</b>       | 적용증: 치기 보철                                      | 54 |
| <b>عربي</b>      | إرشادات الاستخدام: البروتيزات الدنتالية         | 59 |
| <b>Polski</b>    | Instrukcja obsługi: Protetyka Dentystyczna      | 63 |

BioHorizons  
2300 Riverchase Center  
Birmingham, AL 35244 USA  
TOLL FREE 888.246.8338  
TEL 205.967.7880 FAX 205.870.0304  
www.biohorizons.com

Altatec GmbH  
Maybachstrasse 5  
71299 Wimsheim  
Germany  
TEL +49.7044.94450

*The CE mark is valid  
only if it is also printed  
on the product label.*

L02068 Rev C OCT 2022

The symbol table below is for reference only. Refer to product packaging label for applicable symbols.

| Symbol | Symbol Description                                                                                                                                                                                                                                                                     |
|--------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|        | Caution                                                                                                                                                                                                                                                                                |
|        | Electronic Instructions for Use                                                                                                                                                                                                                                                        |
|        | Manufacturer                                                                                                                                                                                                                                                                           |
|        | Intra-Lock products carrying the European Conformity (CE) mark fulfill the requirements of the Medical Device Directive 93/42/EEC as amended by Directive 2007/47/EC or the Medical Devices Regulation 2017/745. The CE mark is valid only if it is also printed on the product label. |
|        | Reference/article number                                                                                                                                                                                                                                                               |
|        | Lot/batch number                                                                                                                                                                                                                                                                       |
|        | Do not re-use                                                                                                                                                                                                                                                                          |
|        | Use-by-date                                                                                                                                                                                                                                                                            |
|        | Sterile by gamma irradiation                                                                                                                                                                                                                                                           |
|        | Date of manufacture                                                                                                                                                                                                                                                                    |
|        | Caution: U.S. Federal law restricts these devices for sale, distribution, and use by, or on the order of, a dentist or physician.                                                                                                                                                      |
|        | EU Authorised Representative                                                                                                                                                                                                                                                           |
|        | Do not use if package is damaged                                                                                                                                                                                                                                                       |
|        | Medical Device                                                                                                                                                                                                                                                                         |
|        | Non-sterile                                                                                                                                                                                                                                                                            |
|        | Do not re-sterilize                                                                                                                                                                                                                                                                    |
|        | Home position                                                                                                                                                                                                                                                                          |

This document supersedes all prior revisions. Original language is English.

This document applies to Intra-Lock Dental Abutments and Prosthetic Components. Every product packaging label contains a description of the enclosed device.

## **DESCRIPTION**

Intra-Lock prosthetic components are intended for the restoration of dental implants within the specific indications of each implant system. The label on each prosthetic package contains important product information including whether the prosthetic is supplied sterile or non-sterile.

## **INDICATIONS FOR USE**

Intra-Lock Prosthetics are intended for the restoration of Intra-Lock and Zimmer® Dental Screw-Vent® and Tapered Screw-Vent® implants.

Intra-Lock Angled Multi-unit Abutments contain a unique LIGHT BLUE abutment screw specifically designed for use in those abutments. Multi-unit prosthetics are intended for use as a support for multiple-implant supported prostheses.

Overdenture Abutments are intended to retain a tissue-supported multiple-unit overdenture in the mandible and maxilla.

Titanium temporary abutments are intended for long term temporary restorations (>30 days); hexed for single-unit and non-hexed for multiple-unit restorations.

PEEK prosthetics are intended for short term temporary restorations (<30 days).

Intra-Lock CAD/CAM Abutments are dental abutments placed onto a dental implant to provide support for dental prosthetic restorations. The abutments include: 1) Titanium abutment blanks with a pre-machined implant connection where the upper portion may be custom-milled in accordance with a patient-specific design using CAD/CAM techniques; and 2) Titanium bases with a pre-machined implant connection upon which a CAD/CAM designed superstructure may be fitted to complete a two-piece dental abutment. The abutments include an abutment screw for fixation to the underlying implant. The abutments may be used for single-unit (single-tooth) or multiple-unit (bridges and bars) restorations and are compatible with Intra-Lock implant systems and Zimmer® Dental Screw-Vent® and Tapered Screw-Vent® dental implants with a yellow (3.5mm) internal hex-connection mating platform diameter.

All digitally designed abutments and/or copings for use with Intra-Lock CAD/CAM Abutments are intended to be sent to an Intra-Lock-validated mill center for manufacture. Intra-Lock abutments designed using CAD/CAM techniques must fulfill the Intra-Lock allowable range of design parameters.

## **CONTRAINDICATIONS**

Intra-Lock prosthetics should not be used in patients who have contraindicated systemic or uncontrolled local diseases such as blood dyscrasias, diabetes, hyperthyroidism, oral infections or malignancies, renal disease, uncontrolled hypertension, liver problems, leukemia, severe vascular heart disease, hepatitis, immunosuppressive disorder, pregnancy, collagen and bone diseases. Relative contraindications may include habits such as tobacco use, alcohol consumption, poor oral hygiene, bruxism, nail biting, pencil biting and improper tongue habits depending on severity. Intra-Lock Multi-unit prosthetics are contraindicated for single tooth replacements.

## **DIRECTIONS FOR USE**

Proper surgical procedures and restorative techniques are the responsibility of the medical professional. Each clinician must evaluate the appropriateness of the procedure used based on personal medical training and experience as applied to the patient case at hand. Intra-Lock strongly recommends completion of dental implant courses and strict adherence to the Instructions For Use that accompany Intra-Lock products.

Healing abutments should be hand tightened (10-15 Ncm).

Abutment screws must be torqued to 30 Ncm to prevent screw loosening. Intra-Lock Multi-unit Abutments, Custom Milled Abutments, and Titanium Base Abutments contain a LIGHT BLUE abutment screw specifically designed for use in those abutments. Use of the LIGHT BLUE abutment screws with incompatible abutments may result in abutment and/or abutment screw failure.

Intra-Lock Multi-Unit Prosthetic Screws must be torqued to 15 Ncm to prevent screw loosening.

Intra-Lock Multi-Unit Prosthetic Component design parameters are as follows:

- Minimum post height should be equal to or greater than 4mm.
- Reduction of wall thickness is not permitted.
- Angle correction is not permitted. Angled Multi-unit Abutments shall be used for angle correction.

Yellow prosthetic connection compatibility with Zimmer® Dental Screw-Vent® and Tapered Screw-Vent® implants is provided in the following table.

| Platform                  | Zimmer Tapered ScrewVent Implants† |                 | Zimmer ScrewVent Implants* |      |
|---------------------------|------------------------------------|-----------------|----------------------------|------|
| Yellow Connection (3.5mm) | TSVBx<br>TSV4Bx<br>TSVT4Bx         | TSVHx<br>TSV4Hx | TSVMBx<br>SVLx<br>SVMHx    | SVHx |

\*Where variable x = implant length.

Secure Custom Cast Abutments to the supplied polishing protector to protect the implant/abutment interface during laboratory procedures.

For CAD/CAM abutments, a custom abutment design file is created by the customer using a 3Shape abutment library. The abutment design file is converted to a stereolithography (.stl) file using 3Shape software. The .stl file is converted to a numerical control (.nc) file using an appropriate software application. After the custom abutment design is uploaded, the milling machine software executes the necessary commands to mill the final abutment. The digitally designed file is to be sent to a registered and listed Intra-Lock contract manufacturer for milling. Intra-Lock CAD/CAM Abutments are compatible with commercially available dental CAD/CAM systems, such as 3Shape.

Design limitation parameters for the Custom Milled Abutments are included in the following table.

| Feature                     | Custom Milled Platforms   |                             |
|-----------------------------|---------------------------|-----------------------------|
|                             | Grey, Yellow              |                             |
|                             | Min                       | Max                         |
| Wall thickness              | 0.3mm                     | 3.7mm                       |
| Gingival Margin Diameter    | 3.3mm, 3.8mm respectively | 9.5mm                       |
| Gingival Margin Height      | 0.3mm, 0.2mm respectively | 12.3mm, 12.5mm respectively |
| Length of the Abutment Post | *                         | 12.3mm, 12.5mm respectively |
| Taper of the Abutment       | 0°                        | 90°                         |
| Angulation                  | 0°                        | 30°                         |

\* Minimum allowable abutment post dimensions are 3.0mm for screw retained abutment and 4.0mm for cement retained abutment.

Titanium base abutments – Intra-Lock recommends the use of Metoxit Z-CAD® dental ceramic (or equivalent) to fabricate the custom zirconia superstructure and RelyX Unicem (or equivalent) to affix the zirconia superstructure to the titanium base.

Design limitation parameters for the custom zirconia superstructure portion of the Titanium Base Abutments are included in the following table.

| Feature                                      | Titanium Base Platforms   |     |
|----------------------------------------------|---------------------------|-----|
|                                              | Grey, Yellow              |     |
|                                              | Min                       | Max |
| <b>Wall Thickness</b>                        | 0.4mm                     | **  |
| <b>Gingival Margin Diameter</b>              | 3.8mm, 3.8mm respectively | **  |
| <b>Length of the Abutment Post</b>           | 4.0mm                     | **  |
| <b>Taper</b>                                 | Variable based on anatomy | **  |
| <b>Angulation of the Final Titanium Base</b> | 0°                        | 20° |

\*\*Feature maximum sizes are based on patient anatomy and clinician discretion.

## WARNINGS AND PRECAUTIONS

Clinician judgment, as related to individual patient presentations, must always supersede recommendations in any Intra-Lock Instructions for Use (IFU). Clinicians are responsible for understanding the appropriate technical use of Intra-Lock prosthetic components. Additional technical information is available upon request from Intra-Lock, or may be viewed and/or downloaded at [www.intra-lock.com](http://www.intra-lock.com). Contact Intra-Lock Customer Care or your local representative with any questions you have regarding specific IFU. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the EU member state in which the clinician and/or patient is established.

Dental implants can break in function for a number of reasons including overloading due to improper occlusion, metal fatigue, and over-tightening of the implant during insertion. Potential causes of abutment fracture include, but are not limited to: casting titanium above 2010°F (1099°C), inadequate implant support when attached to periodontically compromised teeth, non-passive fit of superstructure, overloading due to improper occlusion, incomplete seating of cemented abutments, and excessive cantilevering of pontics.

If any modifications are made to the implant/abutment interface, the abutment may not properly interface with the implant. The FDA considers the manner of the implant/abutment interface a medical device company subject to FDA rules and regulations.

Use of the LIGHT BLUE abutment screws (packaged with Angled Multi-unit Abutments, Custom Milled Abutments, and Titanium Bases) with incompatible components may result in abutment and/or abutment screw failure.

Small diameter implants and angled abutments are not recommended for use in the posterior region of the mouth.

Intra-Lock Multi-unit Abutments require the addition of a coping or cylinder for single and multiple-unit restorations. The coping or cylinder gingival collar, angulation, or wall thickness shall not be modified. The coping or cylinder post height shall not be less than 4mm for single-unit restorations.

The Intra-Lock CAD/CAM Abutments have not been evaluated for safety and compatibility in the MR environment. They have not been tested for heating, migration, or image artifact in the MR environment. The safety of Intra-Lock CAD/CAM Abutments in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

Prosthetics are single patient use only. To eliminate the risk of cross-patient contamination re-use should not be attempted. Intra-Lock assumes no responsibility for attempted re-use or re-sterilization between patients.

Devices should not be used in patients with a known allergy or sensitivity to the device material.

## COMPLICATIONS AND ADVERSE EFFECTS

The risks and complications with prosthetic components and implants include, but are not limited to: (1) allergic reaction(s) to implant and/or abutment material; (2) implant and/or abutment breakage; (3) abutment screw and/ or retaining screw loosening; (4) infection requiring revision of the dental implant; (5) nerve damage that could cause permanent weakness, numbness, or pain; (6) histologic responses possibly involving macrophages and/ or fibroblasts; (7) formation of fat emboli; (8) implant loosening requiring revision surgery; (9) maxillary sinus perforation; (10) labial or lingual plate perforation; and (11) bone loss possibly resulting in revision or removal.

## HANDLING AND STERILIZATION

Always handle the product with powder-free gloves and avoid contact with hard objects that may damage the surface. If the product is supplied sterile, it should be considered sterile unless the package has been opened or damaged. Using accepted sterile technique, remove product from the package only after the correct size has been determined and the surgical site has been prepared. Non-sterile product must be cleaned and sterilized (including CAD/CAM abutments after abutment preparation and before final placement in the restorative site). The following cleaning protocol must be used:

- 1) Prepare a detergent bath in a sterile container using a broad-spectrum cleaning or disinfecting agent such as Hu-Friedy's Enzymax® per the manufacturer's recommendations.
- 2) Brush the product to remove visible debris using a soft bristled brush moistened with the prepared detergent solution.
- 3) Thoroughly rinse product under running utility tap water.
- 4) Place product in the sterile container filled with the prepared detergent solution and sonicate for two (2) minutes minimum.
- 5) Thoroughly rinse product under running utility tap water.
- 6) Spray or wipe product with 70% IPA.
- 7) Blot product dry with clean lint free cloth.

For sterilization of non-sterile product, place product in a FDA cleared sterilization bag or wrap and run through one of the following qualified sterilization cycles:

| Sterilization Cycles           |                            |                             |                            |                            |
|--------------------------------|----------------------------|-----------------------------|----------------------------|----------------------------|
| Reference:                     | ANSI/AAMI TIR12:2010       | ANSI/AAMI TIR12:2010        | ANSI/AAMI TIR12:2010       | UK HTM 01-01 Part C:2016   |
| Type:                          | Gravity Steam              | Gravity Steam               | Prevacuum Steam            | Prevacuum Steam            |
| Exposure Time and Temperature: | 3 minutes at 121°C (250°F) | 15 minutes at 132°C (270°F) | 4 minutes at 132°C (270°F) | 3 minutes at 134°C (273°F) |
| Minimum Dry Time:              | 30 minutes                 | 30 minutes                  | 20 minutes                 | 20 minutes                 |

It is recommended to include a 30-minute cool-down period before removing the product from the sterilization bag or wrap.

Zimmer®, Dental ScrewVent®, and Tapered ScrewVent® are registered trademarks of Zimmer Dental, Inc. Z-CAD® is a registered trademark of Metoxit AG.

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